

LGALS2

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:200 - 1:1000 IHC 1:50 - 1:200

Calculated MW:15kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human LGALS2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HL14;

Catalog #:A2734

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3957

Isotype:IgG

Swiss Prot:P05162

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a soluble beta-galactoside binding lectin. The encoded protein is found as a homodimer and can bind to lymphotoxin-alpha. A single nucleotide polymorphism in an intron of this gene can alter the transcriptional level of the protein, with a resultant increased risk of myocardial infarction.

To place an order, please [Click HERE](#).