

GDF9

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:1000 IHC 1:50 - 1:100

Calculated MW:51kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human GDF9

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GDF9

Catalog #:A2739

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2661

Isotype:IgG

Swiss Prot:O60383

Purity:Affinity purification

For research use only.

Background:

Growth factors synthesized by ovarian somatic cells directly affect oocyte growth and function. Growth differentiation factor-9 (GDF9) is expressed in oocytes and is thought to be required for ovarian folliculogenesis. GDF9 is a member of the transforming growth factor-beta superfamily.

To place an order, please [Click HERE](#).