

HMGB4

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:1000

Calculated MW:22kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human HMGB4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

dJ1007G16.5;

To place an order, please [Click HERE](#).

Catalog #:A2744

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:127540

Isotype:IgG

Swiss Prot:Q8WW32

Purity:Affinity purification

For research use only.