

IGF2BP2

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 66kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human IGF2BP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

IMP2; IMP-2; VICKZ2;

Catalog #: A2749

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10644

Isotype: IgG

Swiss Prot: Q9Y6M1

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the IGF-II mRNA-binding protein (IMP) family. The protein encoded by this gene contains four KH domains and two RRM domains. It functions by binding to the 5' UTR of the insulin-like growth factor 2 (IGF2) mRNA and regulating IGF2 translation. Alternative promoter usage and alternate splicing result in multiple variants encoding different isoforms.

To place an order, please [Click HERE](#).