

MYOZ1

Reactivity:Human Mouse Rat

Tested applications:WB IP

Recommended Dilution:WB 1:500 - 1:2000 IP 1:20 - 1:50

Calculated MW:32kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human MYOZ1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CS-2; FATZ; MYOZ;

Catalog #:A2767

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:58529

Isotype:IgG

Swiss Prot:Q9NP98

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is primarily expressed in the skeletal muscle, and belongs to the myozenin family. Members of this family function as calcineurin-interacting proteins that help tether calcineurin to the sarcomere of cardiac and skeletal muscle. They play an important role in modulation of calcineurin signaling.

To place an order, please [Click HERE](#).