

NME3

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:19kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human NME3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

NDPKC; NDPK-C; NM23H3; DR-nm23; NM23-H3; c371H6.2;

To place an order, please [Click HERE](#).

Catalog #:A2774

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4832

Isotype:IgG

Swiss Prot:Q13232

Purity:Affinity purification

For research use only.