

OASL

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:200 - 1:1000 IHC 1:25 - 1:100

Calculated MW:59kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human OASL

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

OASLd; TRIP14; TRIP-14; p59OASL; p59 OASL; p59-OASL;

To place an order, please [Click HERE](#).

Catalog #:A2775

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8638

Isotype:IgG

Swiss Prot:Q15646

Purity:Affinity purification

For research use only.