

FAM3A

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:200

Calculated MW: 25kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human FAM3A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DLD; 2.19; XAP-7; DXS560S;

Catalog #: A2784

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 60343

Isotype: IgG

Swiss Prot: P98173

Purity: Affinity purification

For research use only.

Background:

This gene encodes a cytokine-like protein. The expression of this gene may be regulated by peroxisome proliferator-activated receptor gamma, and the encoded protein may be involved in the regulation of glucose and lipid metabolism. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).