

RCHY1

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 30kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RCHY1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ZCHY; ARNIP; CHIMP; PIRH2; RNF199; ZNF363; PRO1996;

Catalog #: A2789

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 25898

Isotype: IgG

Swiss Prot: Q96PM5

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene has ubiquitin ligase activity. It mediates E3-dependent ubiquitination and proteasomal degradation of target proteins, including tumor protein 53, histone deacetylase 1, and cyclin-dependent kinase inhibitor 1B, thus regulating their levels and cell cycle progression. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

To place an order, please [Click HERE](#).