

RCHY1

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:1000

Calculated MW:30kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human RCHY1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ZCHY; ARNIP; CHIMP; PIRH2; RNF199; ZNF363; PRO1996;

Background:

The protein encoded by this gene has ubiquitin ligase activity. It mediates E3-dependent ubiquitination and proteasomal degradation of target proteins, including tumor protein 53, histone deacetylase 1, and cyclin-dependent kinase inhibitor 1B, thus regulating their levels and cell cycle progression. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

To place an order, please [Click HERE](#).

Catalog #:A2789

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:25898

Isotype:IgG

Swiss Prot:Q96PM5

Purity:Affinity purification

For research use only.