

PTPN11

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:1000 IHC 1:50 - 1:100

Calculated MW:68kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human PTPN11

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

o

Synonym:

CFC;NS1;SHP2;BPTP3;PTP2C;PTP-1D;SH-PTP2;SH-PTP;

Catalog #:A2793

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5781

Isotype:IgG

Swiss Prot:Q06124

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This PTP contains two tandem Src homology-2 domains, which function as phospho-tyrosine binding domains and mediate the interaction of this PTP with its substrates. This PTP is widely expressed in most tissues and plays a regulatory role in various cell signaling events that are important for a diversity of cell functions, such as mitogenic activation, metabolic control, transcription regulation, and cell migration. Mutations in this gene are a cause of Noonan syndrome as well as acute myeloid leukemia. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).