

RAB17

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:25 - 1:100

Calculated MW:23kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human RAB17

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A2811

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:64284

Isotype:IgG

Swiss Prot:Q9H0T7

Purity:Affinity purification

For research use only.

Background:

The Rab subfamily of small GTPases plays an important role in the regulation of membrane trafficking. RAB17 is an epithelial cell-specific GTPase (Lutcke et al., 1993)

To place an order, please [Click HERE](#).