

RAB18

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:23kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human RAB18

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

WARBM3; RAB18LI1;

Catalog #:A2812

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:22931

Isotype:IgG

Swiss Prot:Q9NP72

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of a family of Ras-related small GTPases that regulate membrane trafficking in organelles and transport vesicles. Knockdown studies in zebrafish suggest that this protein may have a role in eye and brain development. Mutations in this gene are associated with Warburg micro syndrome type 3. Alternatively spliced transcript variants have been found for this gene.

To place an order, please [Click HERE](#).