

RHOB

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:22kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human RHOB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

ab

Synonym:

ARH6; ARHB; RHOH6; MST081; MSTP081;

To place an order, please [Click HERE](#).

Catalog #:A2819

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:388

Isotype:IgG

Swiss Prot:P62745

Purity:Affinity purification

For research use only.