

## RAC1

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**23kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human RAC1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

iqd

**Synonym:**

MIG5; Rac-1; TC-25; p21-Rac1;

**Catalog #:**A2820

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**5879

**Isotype:**IgG

**Swiss Prot:**P63000

**Purity:**Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is a GTPase which belongs to the RAS superfamily of small GTP-binding proteins. Members of this superfamily appear to regulate a diverse array of cellular events, including the control of cell growth, cytoskeletal reorganization, and the activation of protein kinases. Two transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*