

## AVEN

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**39kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic peptide of human AVEN

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

PDCD12;

*To place an order, please [Click HERE](#).*

**Catalog #:**A2843

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**57099

**Isotype:**IgG

**Swiss Prot:**Q9NQS1

**Purity:**Affinity purification

For research use only.