

## CCR7

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**43kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human CCR7

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

BLR2; EBI1; CD197; CDw197; CMKBR7

**Catalog #:**A2857

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**1236

**Isotype:**IgG

**Swiss Prot:**P32248

**Purity:**Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is a member of the G protein-coupled receptor family. This receptor was identified as a gene induced by the Epstein-Barr virus (EBV), and is thought to be a mediator of EBV effects on B lymphocytes. This receptor is expressed in various lymphoid tissues and activates B and T lymphocytes. It has been shown to control the migration of memory T cells to inflamed tissues, as well as stimulate dendritic cell maturation. The chemokine (C-C motif) ligand 19 (CCL19/ECL) has been reported to be a specific ligand of this receptor.

*To place an order, please [Click HERE](#).*