

## CEBPD

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**28kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic peptide of human CEBPD

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

CELF; CRP3; C/EBP-delta; NF-IL6-beta;

**Catalog #:**A2862

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**1052

**Isotype:**IgG

**Swiss Prot:**P49716

**Purity:**Affinity purification

For research use only.

**Background:**

The protein encoded by this intronless gene is a bZIP transcription factor which can bind as a homodimer to certain DNA regulatory regions. It can also form heterodimers with the related protein CEBP-alpha. The encoded protein is important in the regulation of genes involved in immune and inflammatory responses, and may be involved in the regulation of genes associated with activation and/or differentiation of macrophages. The cytogenetic location of this locus has been reported as both 8p11 and 8q11.

*To place an order, please [Click HERE](#).*