

## CHGB

---

**Reactivity:**Human Mouse

**Tested applications:**WB

**Recommended Dilution:**WB 1:100 - 1:500

**Calculated MW:**78kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human CHGB

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

SCG1;

**Catalog #:**A2864

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**1114

**Isotype:**IgG

**Swiss Prot:**P05060

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a tyrosine-sulfated secretory protein abundant in peptidergic endocrine cells and neurons. This protein may serve as a precursor for regulatory peptides.

*To place an order, please [Click HERE](#).*