

CSNK1A1

Reactivity:Human Mouse Rat

Tested applications:WB IF

Recommended Dilution:WB 1:500 - 1:2000 IF 1:20 - 1:100

Calculated MW:39kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human CSNK1A1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CK1; CK1a; CKIa; HLCDDGP1; PRO2975; HEL-S-77p;

To place an order, please [Click HERE](#).

Catalog #:A2868

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1452

Isotype:IgG

Swiss Prot:P48729

Purity:Affinity purification

For research use only.