

KRT16

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 51kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human KRT16

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

K16; PC1; CK16; K1CP; NEPPK; FNEPPK; KRT16A;

Catalog #: A2871

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3868

Isotype: IgG

Swiss Prot: P08779

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the keratin gene family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. Most of the type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains and are clustered in a region of chromosome 17q12-q21. This keratin has been coexpressed with keratin 14 in a number of epithelial tissues, including esophagus, tongue, and hair follicles. Mutations in this gene are associated with type 1 pachyonychia congenita, non-epidermolytic palmoplantar keratoderma and unilateral palmoplantar verrucous nevus.

To place an order, please [Click HERE](#).