

CLEC5A

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:22kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human CLEC5A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MDL1; MDL-1; CLECSF5;

Catalog #:A2873

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:23601

Isotype:IgG

Swiss Prot:Q9NY25

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the C-type lectin/C-type lectin-like domain (CTL/CTLD) superfamily. Members of this family share a common protein fold and have diverse functions, such as cell adhesion, cell-cell signalling, glycoprotein turnover, and roles in inflammation and immune response. The encoded type II transmembrane protein interacts with dnax-activation protein 12 and may play a role in cell activation. Alternative splice variants have been described but their full-length sequence has not been determined.

To place an order, please [Click HERE](#).