

COX18

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:37kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human COX18

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

COX18HS;

Catalog #:A2875

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:285521

Isotype:IgG

Swiss Prot:Q8N8Q8

Purity:Affinity purification

For research use only.

Background:

This gene encodes a cytochrome c oxidase assembly protein. The encoded protein is essential for integral membrane protein insertion into the mitochondrial inner membrane. It is also required for cytochrome c oxidase assembly and activity. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).