

CTGF

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:20 - 1:50

Calculated MW:38kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human CTGF

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

a

Synonym:

CCN2; NOV2; HCS24; IGFBP8;

Catalog #:A2880

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1490

Isotype:IgG

Swiss Prot:P29279

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a mitogen that is secreted by vascular endothelial cells. The encoded protein plays a role in chondrocyte proliferation and differentiation, cell adhesion in many cell types, and is related to platelet-derived growth factor. Certain polymorphisms in this gene have been linked with a higher incidence of systemic sclerosis. [provided by RefSeq, Nov 2009]

To place an order, please [Click HERE](#).