

GJD2

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:2000

Calculated MW:36kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human GJD2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CX36; GJA9;

Catalog #:A2883

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:57369

Isotype:IgG

Swiss Prot:Q9UKL4

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the connexin protein family. Connexins are gap junction proteins which are arranged in groups of 6 around a central pore to form a connexon, a component of the gap junction intercellular channel. The channels formed by this protein allow cationic molecule exchange between human beta cells and may function in the regulation of insulin secretion.

To place an order, please [Click HERE](#).