

AQP10

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:1000

Calculated MW:32kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human AQP10

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AQPA_HUMAN;

Background:

This gene encodes a member of the aquaglyceroporin family of integral membrane proteins. Members of this family function as water-permeable channels in the epithelia of organs that absorb and excrete water. This protein was shown to function as a water-selective channel, and could also permeate neutral solutes such as glycerol and urea.

To place an order, please [Click HERE](#).

Catalog #:A2888

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:89872

Isotype:IgG

Swiss Prot:Q96PS8

Purity:Affinity purification

For research use only.