

CX3CR1

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:44kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human CX3CR1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

V28I;CCRL1;GPR13;CMKDR1;GPRV28;CMKBRL1;

Background:

Fractalkine is a transmembrane protein and chemokine involved in the adhesion and migration of leukocytes. The protein encoded by this gene is a receptor for fractalkine. The encoded protein also is a coreceptor for HIV-1, and some variations in this gene lead to increased susceptibility to HIV-1 infection and rapid progression to AIDS. Four transcript variants encoding two different isoforms have been found for this gene.

To place an order, please [Click HERE](#).

Catalog #:A2890

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1524

Isotype:IgG

Swiss Prot:P49238

Purity:Affinity purification

For research use only.