

## SLC1A6

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:1000

**Calculated MW:**62kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human SLC1A6

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

EAAT4;

*To place an order, please [Click HERE](#).*

**Catalog #:**A2904

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**6511

**Isotype:**IgG

**Swiss Prot:**P48664

**Purity:**Affinity purification

For research use only.