

EGR4

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:25 - 1:100

Calculated MW:62kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human EGR4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NGFIC; NGFI-C; PAT133;

To place an order, please [Click HERE](#).

Catalog #:A2910

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1961

Isotype:IgG

Swiss Prot:Q05215

Purity:Affinity purification

For research use only.