

FOXC1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:100 - 1:1000

Calculated MW: 57kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human FOXC1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ARA; IGDA; IHG1; FKHL7; IRID1; RIEG3; FREAC3; FREAC-3;

Catalog #: A2924

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2296

Isotype: IgG

Swiss Prot: Q12948

Purity: Affinity purification

For research use only.

Background:

This gene belongs to the forkhead family of transcription factors which is characterized by a distinct DNA-binding forkhead domain. The specific function of this gene has not yet been determined; however, it has been shown to play a role in the regulation of embryonic and ocular development. Mutations in this gene cause various glaucoma phenotypes including primary congenital glaucoma, autosomal dominant iridogoniodygenesis anomaly, and Axenfeld-Rieger anomaly.

To place an order, please [Click HERE](#).