

## FOXC1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:100 - 1:1000

**Calculated MW:**57kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human FOXC1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

ARA; IGDA; IHG1; FKHL7; IRID1; RIEG3; FREAC3; FREAC-3;

**Catalog #:**A2924

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**2296

**Isotype:**IgG

**Swiss Prot:**Q12948

**Purity:**Affinity purification

For research use only.

**Background:**

This gene belongs to the forkhead family of transcription factors which is characterized by a distinct DNA-binding forkhead domain. The specific function of this gene has not yet been determined; however, it has been shown to play a role in the regulation of embryonic and ocular development. Mutations in this gene cause various glaucoma phenotypes including primary congenital glaucoma, autosomal dominant iridogoniodysgenesis anomaly, and Axenfeld-Rieger anomaly.

*To place an order, please [Click HERE](#).*