

FOXD1

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:46kDa

Observed MW:Refer to figures

Immunogen:

A synthetic Peptide of human FOXD1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FKHL8; FREAC4; FREAC-4;

Catalog #:A2925

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2297

Isotype:IgG

Swiss Prot:Q16676

Purity:Affinity purification

For research use only.

Background:

This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. Studies of the orthologous mouse protein indicate that it functions in kidney development by promoting nephron progenitor differentiation, and it also functions in the development of the retina and optic chiasm. It may also regulate inflammatory reactions and prevent autoimmunity.

To place an order, please [Click HERE](#).