

## FOXD1

**Reactivity:** Human

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 46kDa

**Observed MW:** Refer to figures

**Immunogen:**

A synthetic Peptide of human FOXD1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

FKHL8; FREAC4; FREAC-4;

**Catalog #:** A2925

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 2297

**Isotype:** IgG

**Swiss Prot:** Q16676

**Purity:** Affinity purification

For research use only.

**Background:**

This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. Studies of the orthologous mouse protein indicate that it functions in kidney development by promoting nephron progenitor differentiation, and it also functions in the development of the retina and optic chiasm. It may also regulate inflammatory reactions and prevent autoimmunity.

*To place an order, please [Click HERE](#).*