

FOXD3

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:48kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human FOXD3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AIS1; HFH2; VAMAS2; Genesis;

Background:

This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. Mutations in this gene cause autoimmune susceptibility 1.

To place an order, please [Click HERE](#).

Catalog #:A2926

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:27022

Isotype:IgG

Swiss Prot:Q9UJU5

Purity:Affinity purification

For research use only.