

FOXL2

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:39kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human FOXL2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BPES; PFRK; POF3; BPES1; PINTO;

Background:

This gene encodes a forkhead transcription factor. The protein contains a fork-head DNA-binding domain and may play a role in ovarian development and function. Mutations in this gene are a cause of blepharophimosis syndrome and premature ovarian failure 3.

To place an order, please [Click HERE](#).

Catalog #:A2932

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:668

Isotype:IgG

Swiss Prot:P58012

Purity:Affinity purification

For research use only.