

## FOXP1

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 75kDa

**Observed MW:** Refer to figures

**Immunogen:**

A synthetic peptide of human FOXP1

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

MFH; QRF1; 12CC4; hFKH1B; HSPC215;

**Catalog #:** A2935

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 27086

**Isotype:** IgG

**Swiss Prot:** Q9H334

**Purity:** Affinity purification

For research use only.

**Background:**

This gene belongs to subfamily P of the forkhead box (FOX) transcription factor family. Forkhead box transcription factors play important roles in the regulation of tissue- and cell type-specific gene transcription during both development and adulthood. Forkhead box P1 protein contains both DNA-binding- and protein-protein binding-domains. This gene may act as a tumor suppressor as it is lost in several tumor types and maps to a chromosomal region (3p14.1) reported to contain a tumor suppressor gene(s). Alternative splicing results in multiple transcript variants encoding different isoforms.

*To place an order, please [Click HERE](#).*