

## GABRR1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**56kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic peptide of human GABRR1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

g

**Catalog #:**A2937

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**2569

**Isotype:**IgG

**Swiss Prot:**P24046

**Purity:**Affinity purification

For research use only.

**Background:**

GABA is the major inhibitory neurotransmitter in the mammalian brain where it acts at GABA receptors, which are ligand-gated chloride channels. GABRR1 is a member of the rho subunit family. Several transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*