

## CCKBR

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**48kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic peptide of human CCKBR

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

GASR; CCK-B; CCK2R;

**Background:**

This gene encodes a G-protein coupled receptor for gastrin and cholecystokinin (CCK), regulatory peptides of the brain and gastrointestinal tract. This protein is a type B gastrin receptor, which has a high affinity for both sulfated and nonsulfated CCK analogs and is found principally in the central nervous system and the gastrointestinal tract. A misspliced transcript variant including an intron has been observed in cells from colorectal and pancreatic tumors.

*To place an order, please [Click HERE](#).*

**Catalog #:**A2941

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**887

**Isotype:**IgG

**Swiss Prot:**P51959

**Purity:**Affinity purification

For research use only.