

## CMTM3

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:1000 IHC 1:50 - 1:100

**Calculated MW:**20kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human CMTM3

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

BNAS2; CKLFSF3

**Catalog #:**A2943

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**123920

**Isotype:**IgG

**Swiss Prot:**Q96MX0

**Purity:**Affinity purification

For research use only.

**Background:**

This gene belongs to the chemokine-like factor gene superfamily, a novel family that is similar to the chemokine and the transmembrane 4 superfamilies of signaling molecules. This gene is one of several chemokine-like factor genes located in a cluster on chromosome 16. Alternatively spliced transcript variants containing different 5' UTRs, but encoding the same protein, have been identified.

*To place an order, please [Click HERE](#).*