

## CMTM2

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**27kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic peptide of human CMTM2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

CKLFSF2;

**Catalog #:**A2944

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**146225

**Isotype:**IgG

**Swiss Prot:**Q8TAZ6

**Purity:**Affinity purification

For research use only.

**Background:**

This gene belongs to the chemokine-like factor gene superfamily, a novel family that links the chemokine and the transmembrane 4 superfamilies of signaling molecules. The protein encoded by this gene may play an important role in testicular development.

*To place an order, please [Click HERE](#).*