

UTS2R

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:42kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human UTS2R

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

UTR; UTR2; GPR14; UR-2-R

Catalog #:A2959

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2837

Isotype:IgG

Swiss Prot:Q9UKP6

Purity:Affinity purification

For research use only.

Background:

This gene codes high affinity receptor for urotensin-2 and urotensin-2B. The activity of this receptor is mediated by a G-protein that activate a phosphatidylinositol-calcium second messenger system

To place an order, please [Click HERE](#).