

HCN1

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 99kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human HCN1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BCNG1; HAC-2; BCNG-1

Catalog #: A2969

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 348980

Isotype: IgG

Swiss Prot: O60741

Purity: Affinity purification

For research use only.

Background:

The membrane protein encoded by this gene is a hyperpolarization-activated cation channel that contributes to the native pacemaker currents in heart and neurons. The encoded protein can homodimerize or heterodimerize with other pore-forming subunits to form a potassium channel.

This channel may act as a receptor for sour tastes.

To place an order, please [Click HERE](#).