

TERT

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:120kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human TERT

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TP2; TRT; EST2; TCS1; hTRT; DKCA2; DKCB4; hEST2; PFBMFT1;

Catalog #:A2979

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7015

Isotype:IgG

Swiss Prot:O14746

Purity:Affinity purification

For research use only.

Background:

Telomerase is a ribonucleoprotein polymerase that maintains telomere ends by addition of the telomere repeat TTAGGG. The enzyme consists of a protein component with reverse transcriptase activity, encoded by this gene, and an RNA component which serves as a template for the telomere repeat. Telomerase expression plays a role in cellular senescence, as it is normally repressed in postnatal somatic cells resulting in progressive shortening of telomeres. Deregulation of telomerase expression in somatic cells may be involved in oncogenesis. Studies in mouse suggest that telomerase also participates in chromosomal repair, since de novo synthesis of telomere repeats may occur at double-stranded breaks.

To place an order, please [Click HERE](#).