

IGFBP1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 28kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human IGFBP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AFBP; IBP1; PP12; IGF-BP25; hIGFBP-1

Catalog #: A2981

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3484

Isotype: IgG

Swiss Prot: P08833

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the insulin-like growth factor binding protein (IGFBP) family and encodes a protein with an IGFBP domain and a thyroglobulin type-I domain. The protein binds both insulin-like growth factors (IGFs) I and II and circulates in the plasma. Binding of this protein prolongs the half-life of the IGFs and alters their interaction with cell surface receptors.

To place an order, please [Click HERE](#).