

KLF5

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:200 - 1:1000 IHC 1:50 - 1:100

Calculated MW:51kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human KLF5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

ak

Synonym:

CKLF; IKLF; BTEB3

Catalog #:A2989

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:688

Isotype:IgG

Swiss Prot:Q13887

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the Kruppel-like factor subfamily of zinc finger proteins. Since the protein localizes to the nucleus and binds the epidermal growth factor response element, it is thought to be a transcription factor.

To place an order, please [Click HERE](#).