

KLF5

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:1000 IHC 1:50 - 1:100

Calculated MW: 51kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human KLF5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

ak

Synonym:

CKLF; IKLF; BTEB3

Catalog #: A2989

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 688

Isotype: IgG

Swiss Prot: Q13887

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the Kruppel-like factor subfamily of zinc finger proteins. Since the protein localizes to the nucleus and binds the epidermal growth factor response element, it is thought to be a transcription factor.

To place an order, please [Click HERE](#).