

LAG3

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:4000

Calculated MW:57kDa

Observed MW:Refer to figures

Immunogen:

Recombinant Protein of human LAG3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CD223;

Catalog #:A2996

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3902

Isotype:IgG

Swiss Prot:P18627

Purity:Affinity purification

For research use only.

Background:

Lymphocyte-activation protein 3 belongs to Ig superfamily and contains 4 extracellular Ig-like domains. The LAG3 gene contains 8 exons. The sequence data, exon/intron organization, and chromosomal localization all indicate a close relationship of LAG3 to CD4.

To place an order, please [Click HERE](#).