

BPIFB1

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:25 - 1:100

Calculated MW:52kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human BPIFB1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LPLUNC1; C20orf114;

Catalog #:A3001

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:92747

Isotype:IgG

Swiss Prot:Q8TDL5

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene may be involved in the innate immune response to bacterial exposure in the mouth, nasal cavities, and lungs. The encoded protein is secreted and is a member of the BPI/LBP/PLUNC protein superfamily. This gene is found with other members of the superfamily in a cluster on chromosome 20.

To place an order, please [Click HERE](#).