

BPIFB2

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:25 - 1:100

Calculated MW:49kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human BPIFB2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RYSR; BPIL1; LPLUNC2; C20orf184; dJ726C3.2;

Background:

This gene encodes a member of the lipid transfer/lipopolysaccharide binding protein (LT/LBP) gene family. It is highly expressed in hypertrophic tonsils. This gene and three other members of the LT/LBP gene family form a cluster on the long arm of chromosome 20.

To place an order, please [Click HERE](#).

Catalog #:A3002

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:80341

Isotype:IgG

Swiss Prot:Q8N4F0

Purity:Affinity purification

For research use only.