

MCMBP

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 73kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MCMBP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MCM-BP; C10orf119;

Catalog #: A3017

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79892

Isotype: IgG

Swiss Prot: Q9BTE3

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein which is a component of the hexameric minichromosome maintenance (MCM) complex which regulates initiation and elongation of DNA. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).